DAVOL INC.

Ventrio™ ST Hernia Patch

Self-Expanding Bioresorbable Coated Permanent Mesh for
Soft Tissue Reconstruction with SorbaFlex™ Memory Technology

Patient Information Leaflet

en	Patient Information Leaflet.....	2
fr	Brochure d'information pour les patients	4
es	Folleto de información para el paciente.....	6
pt	Folheto de informações para o doente.....	8
ru	Листок информации для пациентов	10

DAVOL INC.

Ventrio™ ST Hernia Patch

Self-Expanding Bioresorbable Coated Permanent Mesh for
Soft Tissue Reconstruction with SorbaFlex™ Memory Technology

Patient Information Leaflet

Your doctor will give you a patient implant card for the Ventrio™ ST Hernia Patch. Please keep this card somewhere safe and easy to find. The mesh implant information may be helpful for your future medical care.

INDICATIONS

The Ventrio™ ST Hernia Patch is intended for use in the reinforcement of soft tissue, where weakness exists, in procedures involving soft tissue repair, such as for the repair of hernias.

DEVICE MODELS

5950010, 5950020, 5950030, 5950040, 5950050, 5950060, 5950070, 5950080, 5950090

IMPLANT MATERIAL

Ventrio™ ST Hernia Patch is a self-expanding two-layer mesh, made of polypropylene (PP, < 60%) and Sepramesh™ IP Composite mesh, which are sewn together with polytetrafluoroethylene (PTFE, < 5%) thread, to form a pocket. The Sepramesh™ IP Composite layer is co-knitted using PP and polyglycolic acid (PGA, < 15%) fibers (a common resorbable material). The PGA side is made from a sodium hyaluronate (salt), carboxymethyl cellulose (HA-CMC, < 10%), a common food additive, and polyethylene glycol (PEG, < 15%). PEG is also in toothpaste and skin creams. Between the layers is the SorbaFlex™ Memory Technology and a type of thread that goes away slowly (polydioxanone). Ventrio™ ST Hernia Patch contains D & C Violet N.2 colorant (< 0.1%), which is found in the PDO and PGA.

PURPOSE

Ventrio™ ST Hernia Patch is a sterile (germ-free), partially absorbable implant for the reconstruction of soft tissue defects in adults.

WARNINGS AND PRECAUTIONS

The safety and effectiveness of Ventrio™ ST Hernia Patch has not been evaluated in:

- Pregnant or breastfeeding women
- Infants or children
- Cardiovascular tissue (heart or blood vessels)
- Cancer in the abdomen (belly).

POTENTIAL ADVERSE (UNFAVORABLE) REACTIONS

1. Seroma: An accumulation and pocket of clear fluid. This may result in unscheduled drainage or reoperation.
2. Adhesions: Adhesions form when internal scar tissue joins tissues together that are not usually joined. Adhesions to the bowel (the intestine) may lead to bowel obstruction (blockage). This can cause the need for multiple invasive procedures, prolonged hospitalization, long-term nursing care, intravenous nutrition, and/or staged surgery to fix the bowel.
3. Hematoma: An accumulation of blood.
4. Pain (including chronic or long-term pain): There is a possibility of pain which may affect quality of life and everyday functions, such as driving, walking up the stairs, and/or dressing oneself. Pain may continue and may require a second operation to remove the mesh.
5. Inflammation: Swelling, redness, heat, and/or pain produced in an area of your surgery.
6. Infection: All surgery includes the risk of an infection. You may experience fever, chills, pain, reddening, swelling and increased temperature near the mesh. This may require antibiotics and/or additional diagnostic and invasive procedures, including the chance for reoperation and/or removal of the mesh.
7. Foreign Body reaction: Body reaction to a new material causing swelling, redness, heat, and/or pain around the area of surgery.
8. Fibrosis: Extra tissue formed in the surgical area that can be firm like a scar.
9. Erythema: Redness of the skin or tissue related to infection or inflammation.
10. Mesh migration: Mesh moves from where your doctor placed it during surgery.
11. Mesh erosion: Erosion is the mesh wearing into an organ that may result in the need for multiple invasive procedures, prolonged hospitalization,

- long-term nursing care, and/or intravenous nutrition, and/or staged surgery to address the bowel (or another organ).
12. Mesh extrusion: Mesh poking through the skin.
 13. Fistula formation: A fistula is an abnormal connection or passageway between an organ, vessel, or tissue. Fistula formation may result in the need for multiple invasive procedures, prolonged hospitalization, long-term nursing care, intravenous nutrition, and/or staged surgery to address the bowel (or another organ).
 14. Allergic reaction: An allergic reaction to the mesh may be difficult to diagnose in the postoperative period after surgery due to distracting pain, fever, and the use of sedatives and narcotics. These symptoms, if not recognized as an allergic reaction, may result in additional diagnostic and invasive procedures up to and including the need for reoperation and explanting (removing) the mesh.
 15. Recurrence: Return or reappearance of your hernia or soft tissue defect.
 16. If the memory ring is damaged during mesh placement, it could poke a hole in the bowel or skin and get infected.

DEVICE LIFETIME

Ventrio™ ST Hernia Patch has both absorbable and nonabsorbable components (parts). The absorbable part resorbs in less than 30 days. The PGA is essentially resorbed within 24 weeks, while the permanent polypropylene fibers continue to provide support to the reinforced area. PDO is absorbed within 24-32 weeks.

FOLLOW UP

Please contact your doctor if you have pain, reopening of your wound, swelling, or redness at or around your incision, or an allergic reaction.

For additional Ventrio™ ST Hernia Patch details, please go to bd.com.

Bard, Davol, and Ventrio are trademarks and/or registered trademarks of C. R. Bard, Inc.

Copyright © 2025 C. R. Bard, Inc. All rights reserved.

SYMBOL DEFINITIONS

	Catalogue number		Medical device
	Batch code		Date
	Manufacturer: may contain the Date of Manufacture (YYYY-MM-DD)		Hospital
	Patient identification		Health care center or doctor
	Unique device identifier		Patient information website

 Davol Inc.
 Subsidiary of C. R. Bard, Inc.
 100 Crossings Boulevard
 Warwick, RI 02886 USA
 1-401-825-8300 • 1-800-556-6275
bd.com

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Copyright © 2025 C. R. Bard, Inc. Reservados todos los derechos.
Copyright © 2025 C. R. Bard, Inc. Todos os direitos reservados.